

Product Datasheet

CLK1 antibody (orb131749)

Description	Rabbit polyclonal to CLK1. Dual specificity protein kinase CLK1 belongs to the protein kinase superf...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human CLK1
Preservatives	0.01% NaN ₃
Concentration	100 µg/ml BSA
Storage	-20°C
Note	For research use only
Application notes	WB: The apparent protein size on WB may be different from the calculated M.W. due to modifications. Experiment Notes: Rabbit polyclonal antibodies were produced by immunizing animals with a GST-fusion protein containing the N-terminal region of human CLK1.
Clonality	Polyclonal
Purity	Affinity purification
MW	57
Uniprot ID	P49759
Dilution Range	WB:1:1000-1:3000
Expiration Date	12 months from date of receipt.