


Product Datasheet

NUDT10 antibody (orb1317517)

Description	Mouse monoclonal antibody to NUDT10
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human NUDT10 (NP_694853) produced in E.coli.
Target	NUDT10
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	18.3 kDa
Uniprot ID	Q8NFP7
NCBI	NM_153183
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...