

Product Datasheet

TET3 antibody (orb1317586)

Description	Mouse monoclonal antibody to TET3
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 241-568 of human TET3 (NP_659430) produced in E.coli.
Target	TET3
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	O43151
NCBI	NM_144993
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...