

Product Datasheet

GLI1 antibody (orb1317619)

Description	Mouse monoclonal antibody to GLI1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 755-1106 of human GLI1 (NP_005260) produced in E.coli.
Target	GLI1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	117.7 kDa
Uniprot ID	P08151
NCBI	NM_005269
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...