


Product Datasheet

MDA5 (IFIH1) antibody (orb1317718)

Description	Mouse monoclonal antibody to MDA5 (IFIH1)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human IFIH1 (NP_071451) produced in HEK293T cell.
Target	IFIH1
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	116.5 kDa
Uniprot ID	Q9BYX4
NCBI	NM_022168
Dilution Range	WB 1:500~2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...