


Product Datasheet

ZNF34 antibody (orb1317748)

Description	Mouse monoclonal antibody to ZNF34
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 117-405 of human ZNF34 (NP_085057) produced in E.coli.
Target	ZNF34
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	63.9 kDa
Uniprot ID	Q8IZ26
NCBI	NM_030580
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...