

Product Datasheet

SRSF7 antibody (orb131776)

Description	Rabbit polyclonal to SRSF7. Serine/arginine-rich splicing factor 7 (SRSF7) belongs to the splicing f...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human SRSF7
Preservatives	0.01% NaN ₃
Concentration	100 µg/ml BSA
Storage	-20°C
Note	For research use only
Application notes	WB: The apparent protein size on WB may be different from the calculated M.W. due to modifications. Experiment Notes: Rabbit polyclonal antibodies were produced by immunizing animals with a GST-fusion protein containing the C-terminal region of human SRSF7.
Clonality	Polyclonal
Purity	Affinity purification
MW	27
Uniprot ID	Q16629
Dilution Range	WB:1:1000-1:3000
Expiration Date	12 months from date of receipt.