


Product Datasheet

MSMB antibody (orb1317787)

Description	Mouse monoclonal antibody to MSMB
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Purified prostate specific antigen (MSMB) from human seminal fluid
Target	MSMB
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	10.7 kDa
Uniprot ID	P08118
NCBI	NM_002443
Dilution Range	WB 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...