


Product Datasheet

CD63 antibody (orb1317802)

Description	Mouse monoclonal antibody to CD63
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC
Immunogen	Full length human recombinant protein of human CD63 (NP_001771) produced in HEK293T cell.
Target	CD63
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	25.5 kDa
Uniprot ID	P08962
NCBI	NM_001780
Dilution Range	IHC 1:150, FLOW 1:50
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...