


Product Datasheet

Von Willebrand Factor (VWF) antibody (orb1317838)

Description	Mouse monoclonal antibody to Von Willebrand Factor (VWF)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 764-1263 of human VWF (NP_000543) produced in SF9 cell.
Target	VWF
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	309.3 kDa
Uniprot ID	P04275
NCBI	NM_000552
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...