

Product Datasheet

TET2 antibody pairs (orb131789)

Description	Rabbit polyclonal to TET2. TET2 antibodies represent a pair of fully characterized antibodies that r...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ChIP, IP, WB
Immunogen	C-terminal region of human TET2 or the N-terminal region of mouse TET2
Preservatives	0.01% NaN ₃
Concentration	100 µg/ml BSA
Storage	-20°C
Note	For research use only
Application notes	WB: The apparent protein size on WB may be different from the calculated M.W. due to modifications. Experiment Notes: Rabbit polyclonal antibodies were produced by immunizing animals with GST-fusion proteins containing either the C-terminal region of human TET2 or the N-terminal region of mouse.
Clonality	Polyclonal
Purity	Affinity purification
MW	224
Uniprot ID	Q6N021, Q4JK59
Dilution Range	WB:1:1000-1:3000
Expiration Date	12 months from date of receipt.