

Product Datasheet

CALCOCO2 antibody (orb131793)

Description	Rabbit polyclonal to CALCOCO2. Calcium-binding and coiled-coil domain-containing protein 2 (CALCOCO2...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human CALCOCO2
Preservatives	0.01% NaN ₃
Concentration	100 µg/ml BSA
Storage	-20°C
Note	For research use only
Application notes	WB: The apparent protein size on WB may be different from the calculated M.W. due to modifications. Experiment Notes: Rabbit polyclonal antibodies were produced by immunizing animals with a GST-fusion protein containing the N-terminal region of human CALCOCO2.
Clonality	Polyclonal
Purity	Affinity purification
MW	52
Uniprot ID	Q13137
Dilution Range	WB:1:1000-1:3000
Expiration Date	12 months from date of receipt.