


Product Datasheet

PLOD2 antibody (orb1317954)

Description	Mouse monoclonal antibody to PLOD2
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 359-640 of human PLOD2 (NP_891988) produced in E.coli.
Target	PLOD2
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	O00469
NCBI	NM_182943
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...