


Product Datasheet

PRRX1 antibody (orb1317978)

Description	Mouse monoclonal antibody to PRRX1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human PRRX1 (NP_073207) produced in E.coli.
Target	PRRX1
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	27.1 kDa
Uniprot ID	P54821
NCBI	NM_022716
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...