


Product Datasheet

RMND5A antibody (orb1317990)

Description	Mouse monoclonal antibody to RMND5A
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-249 of human RMND5A (NP_073617) produced in E.coli.
Target	RMND5A
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	43.8 kDa
Uniprot ID	Q9H871
NCBI	NM_022780
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...