

Product Datasheet

SDAD1 antibody (orb131800)

Description	Rabbit polyclonal to SDAD1. Protein SDA1 homolog (SDAD1) belongs to the SDA1 family and is required ...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human SDAD1
Preservatives	0.01% NaN ₃
Concentration	100 µg/ml BSA
Storage	-20°C
Note	For research use only
Application notes	WB: The apparent protein size on WB may be different from the calculated M.W. due to modifications. Experiment Notes: Rabbit polyclonal antibodies were produced by immunizing animals with a GST-fusion protein containing the N-terminal region of human SDAD1.
Clonality	Polyclonal
Purity	Affinity purification
MW	80
Uniprot ID	Q9NVU7
Dilution Range	WB:1:1000-1:3000
Expiration Date	12 months from date of receipt.