


Product Datasheet

ARL3 antibody (orb1318021)

Description	Mouse monoclonal antibody to ARL3
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 2-182 of human ARL3 (NP_004302) produced in E.coli.
Target	ARL3
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	P36405
NCBI	NM_004311
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...