


Product Datasheet

EHHADH antibody (orb1318027)

Description	Mouse monoclonal antibody to EHHADH
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 496-723 of human EHHADH (NP_001957) produced in E.coli.
Target	EHHADH
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	79.3 kDa
Uniprot ID	Q08426
NCBI	NM_001966
Dilution Range	WB 1:500, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...