


Product Datasheet

DAP3 antibody (orb1318095)

Description	Mouse monoclonal antibody to DAP3
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 161-398 of human DAP3 (NP_387506) produced in E.coli.
Target	DAP3
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	45.4 kDa
Uniprot ID	P51398
NCBI	NM_033657
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...