


Product Datasheet

MRPS11 antibody (orb1318124)

Description	Mouse monoclonal antibody to MRPS11
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human MRPS11 (NP_073750) produced in E.coli.
Target	MRPS11
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	20.4 kDa
Uniprot ID	P82912
NCBI	NM_022839
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...