

Product Datasheet

SAFB antibody (orb131814)

Description	Rabbit polyclonal to SAFB. Scaffold attachment factor B (SAFB) is a DNA-binding protein which has hi...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human SAFB
Preservatives	0.01% NaN ₃
Concentration	100 µg/ml BSA
Storage	-20°C
Note	For research use only
Application notes	WB: The apparent protein size on WB may be different from the calculated M.W. due to modifications. Experiment Notes: Rabbit polyclonal antibodies were produced by immunizing animals with a GST-fusion protein containing the C-terminal region of human SAFB.
Clonality	Polyclonal
Purity	Affinity purification
MW	103
Uniprot ID	Q15424
Dilution Range	WB:1:1000-1:3000
Expiration Date	12 months from date of receipt.