


Product Datasheet

CSNK1G2 antibody (orb1318148)

Description	Mouse monoclonal antibody to CSNK1G2
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 163-415 of human CSNK1G2 (NP_001310) produced in E.coli.
Target	CSNK1G2
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	47.3 kDa
Uniprot ID	P78368
NCBI	NM_001319
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...