

Product Datasheet

GRIN1 antibody (orb131815)

Description Rabbit polyclonal to GRIN1. Glutamate receptor ionotropic, NMDA 1 (GRIN1) belongs to the glutamate-g...

Species/Host Rabbit

Reactivity Human, Mouse

Conjugation Unconjugated

Tested Applications WB

Immunogen C-terminal region of human GRIN1

Preservatives 0.01% NaN₃

Concentration 100 µg/ml BSA

Storage -20°C

Note For research use only

Application notes **WB:** The apparent protein size on WB may be different from the calculated M.W. due to modifications.

Experiment Notes: Rabbit polyclonal antibodies were produced by immunizing animals with a GST-fusion protein containing the C-terminal region of human GRIN1.

Clonality Polyclonal

Purity Affinity purification

MW 105

Uniprot ID **Q05586**

Dilution Range WB:1:1000-1:3000

Expiration Date 12 months from date of receipt.