


Product Datasheet

TIF1 alpha (TRIM24) antibody (orb1318168)

Description	Mouse monoclonal antibody to TIF1 alpha (TRIM24)
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 706-1016 of human TRIM24 (NP_003843) produced in E.coli.
Target	TRIM24
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	O15164
NCBI	NM_003852
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of E.coli lysate (...)