

Product Datasheet

PAK3 antibody (orb131818)

Description	Rabbit polyclonal to PAK3. Serine/threonine-protein kinase PAK 3 (PAK3) belongs to the protein kinas...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of human PAK3
Preservatives	0.01% NaN ₃
Concentration	100 µg/ml BSA
Storage	-20°C
Note	For research use only
Application notes	WB: The apparent protein size on WB may be different from the calculated M.W. due to modifications. Experiment Notes: Rabbit polyclonal antibodies were produced by immunizing animals with a GST-fusion protein containing the middle region of human PAK3.
Clonality	Polyclonal
Purity	Affinity purification
MW	62
Uniprot ID	075914
Dilution Range	WB:1:1000-1:3000
Expiration Date	12 months from date of receipt.