

Product Datasheet

SETDB2 antibody (orb131820)

Description Rabbit polyclonal to SETDB2. Histone-lysine N-methyltransferase SETDB2 belongs to the histone-lysine...

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications WB

Immunogen N-terminal region of human SETDB2

Preservatives 0.01% NaN₃

Concentration 100 µg/ml BSA

Storage -20°C

Note For research use only

Application notes **WB:** The apparent protein size on WB may be different from the calculated M.W. due to modifications.

Experiment Notes: Rabbit polyclonal antibodies were produced by immunizing animals with a GST-fusion protein containing the N-terminal region of human SETDB2.

Clonality Polyclonal

Purity Affinity purification

MW 82

Uniprot ID **Q96T68**

Dilution Range WB:1:1000-1:3000

Expiration Date 12 months from date of receipt.