


Product Datasheet

DOT1L antibody (orb1318206)

Description	Mouse monoclonal antibody to DOT1L
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1008-1424 of human DOT1L (NP_115871) produced in E.coli.
Target	DOT1L
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	164.7 kDa
Uniprot ID	Q8TEK3
NCBI	NM_032482
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...