


Product Datasheet

Progesterone Receptor (PGR) antibody (orb1318219)

Description	Mouse monoclonal antibody to Progesterone Receptor (PGR)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-298 of human PGR (NP_000917) produced in E.coli.
Target	PGR
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	98.8 kDa
Uniprot ID	P06401
NCBI	NM_000926
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...