


Product Datasheet

CDX2 antibody (orb1318248)

Description	Mouse monoclonal antibody to CDX2
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human CDX2 (NP_001256) produced in E.coli.
Target	CDX2
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	33.3 kDa
Uniprot ID	Q99626
NCBI	NM_001265
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...