


Product Datasheet

PARN antibody (orb1318260)

Description	Mouse monoclonal antibody to PARN
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-266 of human PARN (NP_002573) produced in E.coli.
Target	PARN
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	O95453
NCBI	NM_002582
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...