

Product Datasheet

STAT3 antibody (orb131830)

Description	Rabbit polyclonal to STAT3. Signal transducer and activator of transcription 3 (STAT3) is a member o...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human STAT3
Preservatives	0.01% NaN ₃
Concentration	100 µg/ml BSA
Storage	-20°C
Note	For research use only
Application notes	<p>WB: The apparent protein size on WB may be different from the calculated M.W. due to modifications.</p> <p>Experiment Notes: Rabbit polyclonal antibodies were produced by immunizing animals with a GST-fusion protein containing the C-terminal region of human STAT3.</p>
Clonality	Polyclonal
Purity	Affinity purification
MW	88
Uniprot ID	P40763
Dilution Range	WB:1:1000-1:3000
Expiration Date	12 months from date of receipt.