


Product Datasheet

Topoisomerase II alpha (TOP2A) antibody (orb1318305)

Description	Mouse monoclonal antibody to Topoisomerase II alpha (TOP2A)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1100-1531 of human TOP2A(NP_001058) produced in E.coli.
Target	TOP2A
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	174.2 kDa
Uniprot ID	P11388
NCBI	NM_001067
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...