

Product Datasheet

CAPN1 antibody (orb131835)

Description Rabbit polyclonal to CAPN1. Calpain-1 catalytic subunit (CAPN1) belongs to the peptidase C2 family a...

Species/Host Rabbit

Reactivity Human, Mouse

Conjugation Unconjugated

Tested Applications WB

Immunogen C-terminal region of human CAPN1

Preservatives 0.01% NaN₃

Concentration 100 µg/ml BSA

Storage -20°C

Note For research use only

Application notes **WB:** The apparent protein size on WB may be different from the calculated M.W. due to modifications.

Experiment Notes: Rabbit polyclonal antibodies were produced by immunizing animals with a GST-fusion protein containing the C-terminal region of human CAPN1.

Clonality Polyclonal

Purity Affinity purification

MW 82

Uniprot ID **P07384**

Dilution Range WB:1:1000-1:3000

Expiration Date 12 months from date of receipt.