

Product Datasheet

PAK7 antibody (orb131838)

Description	Rabbit polyclonal to PAK7. Serine/threonine-protein kinase PAK 7 belongs to the protein kinase super...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of human PAK7
Preservatives	0.01% NaN ₃
Concentration	100 µg/ml BSA
Storage	-20°C
Note	For research use only
Application notes	WB: The apparent protein size on WB may be different from the calculated M.W. due to modifications. Experiment Notes: Rabbit polyclonal antibodies were produced by immunizing animals with a GST-fusion protein containing the middle region of human PAK7.
Clonality	Polyclonal
Purity	Affinity purification
MW	81
Uniprot ID	Q9P286
Dilution Range	WB:1:1000-1:3000
Expiration Date	12 months from date of receipt.