


Product Datasheet

AFAP (AFAP1) antibody (orb1318392)

Description	Mouse monoclonal antibody to AFAP (AFAP1)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 100-386 of human AFAP1 (NP_940997) produced in E.coli.
Target	AFAP1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	80.5 kDa
Uniprot ID	Q8N556
NCBI	NM_198595
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...