

Product Datasheet

5-OxoETE peptide (orb13184)

Description	5-OxoETE peptide
Conjugation	Unconjugated
Tested Applications	HPLC
Target	5-OxoETE
Form/Appearance	Lyophilized
Storage	Upon reconstitution, aliquot and store at -20°C or -80°C. for long term storage, Avoid repeated freeze/thaw cycles.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb385608 (5-OxoETE antibody)
Expiration Date	6 months from date of receipt.