


Product Datasheet

CD4 antibody (orb1318419)

Description	Mouse monoclonal antibody to CD4
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Full length human recombinant protein of human CD4 (NP_000607) produced in HEK293T cell.
Target	CD4
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	48.3 kDa
Uniprot ID	P01730
NCBI	NM_000616
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...