


Product Datasheet

OXSM antibody (orb1318434)

Description	Mouse monoclonal antibody to OXSM
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 78-343 of human OXSM (NP_060367) produced in E.coli.
Target	OXSM
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	48.7 kDa
Uniprot ID	Q9NWU1
NCBI	NM_017897
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...