

Product Datasheet

GOLGA5 antibody (orb131845)

Description Rabbit polyclonal to GOLGA5. Golgin subfamily A member 5 (GOLGA5) is involved in maintaining Golgi s...

Species/Host Rabbit

Reactivity Human, Mouse

Conjugation Unconjugated

Tested Applications WB

Immunogen C-terminal region of human GOLGA5

Preservatives 0.01% NaN₃

Concentration 100 µg/ml BSA

Storage -20°C

Note For research use only

Application notes **WB:** The apparent protein size on WB may be different from the calculated M.W. due to modifications.

Experiment Notes: Rabbit polyclonal antibodies were produced by immunizing animals with a GST-fusion protein containing the C-terminal region of human GOLGA5.

Clonality Polyclonal

Purity Affinity purification

MW 83

Uniprot ID **Q8TBA6**

Dilution Range WB:1:1000-1:3000

Expiration Date 12 months from date of receipt.