


## Product Datasheet

### VILIP1 (VSNL1) antibody (orb1318527)

|                            |                                                                                                                      |
|----------------------------|----------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to VILIP1 (VSNL1)                                                                          |
| <b>Species/Host</b>        | Mouse                                                                                                                |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                                    |
| <b>Conjugation</b>         | Unconjugated                                                                                                         |
| <b>Tested Applications</b> | WB                                                                                                                   |
| <b>Immunogen</b>           | Human recombinant protein fragment corresponding to amino acids 2-191 of human VSNL1 (NP_003376) produced in E.coli. |
| <b>Target</b>              | VSNL1                                                                                                                |
| <b>Note</b>                | For research use only                                                                                                |
| <b>Isotype</b>             | IgG1                                                                                                                 |
| <b>Clonality</b>           | Monoclonal                                                                                                           |
| <b>MW</b>                  | 22 kDa                                                                                                               |
| <b>Uniprot ID</b>          | <b>P62760</b>                                                                                                        |
| <b>NCBI</b>                | <b>NM_003385</b>                                                                                                     |
| <b>Dilution Range</b>      | WB 1:2000                                                                                                            |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                      |



Western blot analysis of HEK293T cells w...