


Product Datasheet

Hamartin (TSC1) antibody (orb1318545)

Description	Mouse monoclonal antibody to Hamartin (TSC1)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 836-1164 of human TSC1 (NP_000359) produced in E.coli.
Target	TSC1
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	129.6 kDa
Uniprot ID	Q92574
NCBI	NM_000368
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...