


Product Datasheet

HMG2L1 (HMGXB4) antibody (orb1318553)

Description	Mouse monoclonal antibody to HMG2L1 (HMGXB4)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-348 of human HMGXB4 (NP_005478) produced in E.coli.
Target	HMGXB4
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	Q9UGU5
NCBI	NM_005487
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...