


Product Datasheet

WNT3 antibody (orb1318590)

Description	Mouse monoclonal antibody to WNT3
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 128-355 of human WNT3 (NP_110380) produced in E.coli.
Target	WNT3
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P56703
NCBI	NM_030753
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...