


Product Datasheet

FOXP1 antibody (orb1318603)

Description	Mouse monoclonal antibody to FOXP1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC
Immunogen	Full length human recombinant protein of human FOXP1 (NP_116071) produced in SF9 cell.
Target	FOXP1
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	75.1 kDa
Uniprot ID	Q9H334
NCBI	NM_032682
Dilution Range	IHC 1:150, FLOW 1:50
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...