


Product Datasheet

PIK3CD antibody (orb1318610)

Description	Mouse monoclonal antibody to PIK3CD
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 286-610 of human PIK3CD (NP_005017) produced in E.coli.
Target	PIK3CD
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	119.3 kDa
Uniprot ID	O00329
NCBI	NM_005026
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of E.coli lysate (...)