


Product Datasheet

Kallikrein 8 (KLK8) antibody (orb1318638)

Description	Mouse monoclonal antibody to Kallikrein 8 (KLK8)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 33-260 of human KLK8 (NP_009127) produced in SF9 cell.
Target	KLK8
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	O60259
NCBI	NM_007196
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...