

Product Datasheet

PI 3 Kinase Class 2A (PIK3C2A) antibody (orb1318641)

Description Mouse monoclonal antibody to PI 3 Kinase Class 2A (PIK3C2A)

Species/Host Mouse

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications WB

Immunogen Human recombinant protein fragment corresponding to amino acids 230-560 of human PIK3C2A (NP_002636) produced in E.coli.

Target PIK3C2A

Note For research use only

Isotype IgG2b

Clonality Monoclonal

MW 190.5 kDa

Uniprot ID [O00443](#)

NCBI [NM_002645](#)

Dilution Range WB 1:2000

Expiration Date 12 months from date of receipt.

