


Product Datasheet

KIBRA (WWC1) antibody (orb1318700)

Description	Mouse monoclonal antibody to KIBRA (WWC1)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human WWC1 (NP_056053) produced in HEK293T cell.
Target	WWC1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	125.1 kDa
Uniprot ID	Q8IX03
NCBI	NM_015238
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...