


Product Datasheet

Thymidylate Synthase (TYMS) antibody (orb1318705)

Description	Mouse monoclonal antibody to Thymidylate Synthase (TYMS)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human TYMS (NP_001062) produced in E.coli.
Target	TYMS
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	35.5 kDa
Uniprot ID	P04818
NCBI	NM_001071
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...