


Product Datasheet

CD23 (FCER2) antibody (orb1318706)

Description	Mouse monoclonal antibody to CD23 (FCER2)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 48-321 of human FCER2 (NP_001993) produced in Insect.
Target	FCER2
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	36.3 kDa
Uniprot ID	P06734
NCBI	NM_002002
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...